



POCKET SLIT LAMP

PORTABLE | VISION TEST | OPHTHALMIC CHECK | HIGH CLEAR



Specification:

Magnification: 5X
 Slit Length: 5mm-11mm
 Min. Slit With: 0.1mm
 Working Distance: 12-32mm
 Power Supply: 2XAA Battery
 Light Source: 3V/1W LED Lamp
 Battery Life: More Than 3 Hours
 Net Weight: 120g

Model	Discription	Photo
S1501	Round Spot With White Lamp	●
S1502	Slit Light With White Lamp	
S1503	Round Spot With Cobalt Bule Filter	●
S1504	Slit Light With Cobalt Bule Filter	
S1505	Round Spot With Red Free Filter	●
S1506	Slit Light With Red Free Filter	