



TABLE TOP STEAM STERILIZER

Characteristics:

1. **Sterilizing course: Automatic sterilization controlled by computer, easy to operate.**
2. Maximum temperature: up to 134 °C, suit for 4-6 minutes rapidly sterilizing.
3. Safe protection of water lacking.
4. **With over-temperature, over-pressure auto-protect device**
5. After sterilization it will automatically cut the power and then alarm.
6. Fully stainless steel structure.

Technical data:

MODEL TECHNICAL DATA	CCS-T24J (YXQ. DY. 250A-I 25)
Sterilizing chamber volume	24L (ϕ 270×410 mm)
working pressure	0.22Mpa
working temperature	134℃
Adjustment of temperature	105-134℃
Timer	0-60 mi n
Heat average	$\leq \pm 1^{\circ}\text{C}$
Power	1.5 KW / AC220V 50Hz
Sterilizing plate	340 × 200 × 30 mm (3pi ece)
Dimension	560 × 350 × 500 mm
Package dimension	620 × 410 × 570 mm
G.W/N. W	20/18 kg