

C2601016015

Full-digital Ultrasound

Full-digital Ultrasound
Diagnostic System expert



Multifunctional



Language

Built in 8 languages, such as Chinese, English, French, Spanish, etc.



Department

Built in 8 languages, such as Chinese, English, French, Spanish, etc.



Race

Formula for measuring different race. Select standard obstetric algorithms based on local ethnicity.



Report page

One click report generation page. Automatically import diagnostic images and measurement data.

Convenient and vivid

ARM chip architecture, with more stable and powerful functions.

32 channel signal processing for richer and clearer images.

15 inch 1024 * 768 resolution high-definition LCD screen.

Dual probe activation socket, VGA interface, USB2.0 interface, VIDEO interface.

Built in large capacity lithium battery with a range of 2.5 to 5 hours.



2.

Convenient and vivid



ARM chip architecture, with more stable and powerful functions.

32 channel signal processing for richer and clearer images.

15inch 1024 * 768 resolution high-definition LCD screen.

Dual probe activation socket,

VGA interface, USB2.0 interface, VIDEO interface.

Built in **lithium battery** with a range of **2.5 to 5 hours.**

3.

Body marks

97 Body marks




Formula for measuring different race



5. Applications

Formula for measuring different races



4.

Applications



Formula for measuring different races



Ulna	EFW	GS	< Obstetrics
Radius	Last menstruation	BPD	Distance
Tibia	Amniotic Fluid	HC	Circle - Ellipse
Fibula	Amniotic Fluid Index	AC	Circle - Trajectory
Clavicle	Placenta	FL	Area - Ellipse
Vertebra	Th C	CRL	Area - Trajectory
Middle phalanx	Umbilical vein	Yolk Sac	Volume
Foot length	FK	TTD	Angle
Ear length	Umbilical Cord	AFTD	Distance Ratio
Trunk	Cervix length	CTD	Area Ratio
Chest diameter	Neck fold	NT	
Fetal Heart circum	Fetal orbit	OFD	
Fetal heart rate	OOD	Cist mag	
	IOD	LVW	
	Humerus	Falx Cerebri	


6.

PW FUNCTION



7.

Different Applications





8.

Different Applications



9.


Pesudeo Color



10. Editable Report Page



KEY WORDS







Ultrasonic parts : Superficial Ultrasound No :

Name : Gender: Female Age :

Registered No : Applying physician :

Ultrasound image :



Measurement result & Describe :

Right THY L :00 .90 cm , Right THY H :00 .97 cm , Right THY W :01 .37 cm ;

ffffffffffbnm bvf dv

Summary :

csdsbnbnbn

Diagnostician : Reviewer :

For reference of clinicians only ! Audit date :

11. Strong Packing

**Alloy case
PACKED**

**Carton Size:40cm*21cm*42cm
GW:7.9kgs NW:3.5KGS
CBM:0.036**

