



- Ideal for microfiltration and quick spin-downs
- Easy-to-use with dual start/stop function
- 7000rpm fixed speed(CCL-D1008)
- 5000rpm fixed speed(CCL-D1008E)
- Compatible with 2mLx8 and PCR8x2 rotors
- Quiet running, noise level ≤ 45 dB
- Rotor change without tools
- Accelerates and brakes in seconds
- Available in two colors i.e., green, blue

Specifications	CCL-D1008
Max. Speed	7000rpm
Max. RCF	2680 \times g
Rotor Capacity	0.2/0.5/1.5/2.0mL \times 8; 0.2mL \times 16 PCR strips or 0.2mL \times 2 PCR 8 strips
Run Time	Continuous operation
Driving Motor	DC motor
Power	AC110-240V/50Hz/60Hz 20W
Noise Level	≤ 45 dB
Dimension [D \times W \times H]	160x170x122mm

Weight	0.5kg
Certification	CE cTÜVus FCC MCA