



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26 ,Tianhe,Guangzhou, China

Full Automatic Uncap Centrifuge Model: CCL-Z4WS



Features:

1. Table-type uncapped centrifuges, which adopt full steel structure and stainless chamber with small size and low noise
 2. Solved the problem of blood re-mixed by hand-uncovered, which greatly improved work efficiency and completely eradicated bacteria infection and prevent reagents from human body.
 3. It can be done once of the multi-functional of uncap and separation.
- It is the replacement product for vacuum tube centrifuge.

Technical parameters

Model	CCL-Z4WS		
Max Speed	4000r/min		
Rotor	Capacity	Max Speed	Capacity
Size	Small	Medium	Large
Capacity	32×5/7ml	48×5/7ml 72×5/7	96×5/7ml
RCF	2800×g	3050×g	3350×g
Power	450W	750W	750W
Motor	AC Frequency Motor, Microprocessor Control		
Timer Range	0~99min		
Noise	≤65dB(A)		
Power Supply	AC220V&110V 50Hz 5A		
Dimension	580×450×370mm(L×W×H)		



Skype: chinacaremed

Whatsapp: +86 139 2648 0041

Facebook: China care medical

Email: joseph@chinacaremedical.com

Wechat: joseph20130708

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguang South Road No 26, Tianhe, Guangzhou, China

Weight

39kg

52kg

100kg

Rotors

No.1 rotor

Capacity: 32×5/7/10ml, Max speed: 4000 r/m, Max RCF:2800xg

No.2 rotor

Capacity: 48×5/7/10ml, Max speed: 4000 r/m, Max RCF:2800xg

No.3 rotor

Capacity: 72×5/7/10ml, Max speed: 4000 r/m, Max RCF:3050xg

No.4 rotor

Capacity: 96×5/7/10ml, Max speed: 4000 r/m, Max RCF:3680xg