



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Dagan South Road No 26 ,Tianhe,Guangzhou, China

Portable Biphasic AED Defibrillator

Portable Defibrillator, Defibrillator, AED Defibrillator

CCD-AED7000



Specifications:

- 1.three-step defibrillation process
- 2.two-button operation
- 3.extensive voice and visual prompts for the operator
- 4.biphasic energy output
- 5.lock-out protection to prevent inadvertent defibrillation

Battery:

Category	Specification
Part No	CR123A-4x2
Non-rechargeable	12V d.c. 2.8Ah
Capacity	100 discharges at 200 Joules or 120 discharges at 150 Joules
Shelf Life (25°C±15°C)	5 years

Physical:

Category	Specification
Dimensions	303 x 216 x 89 mm
Weight	2.0 kg
Operating Temperature	0°C to 40°C
Operating Humidity	Relative humidity between 30% and 95% (non-condensing)
Storage Temperature(without battery)	-20°C to 55°C
Storage Humidity(without battery)	Up to 93% (non-condensing)

Electrical Isolation:



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Dagan South Road No 26, Tianhe, Guangzhou, China

Category	Specification
Power	Unit operates on internal battery only
External Electrical	No external devices are attached to the unit
Risk Current Category	Internally powered equipment with defibrillator-proof BF type patient applied part (as per definition of IEC 60601-1 standard)