



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

Biphasic Defi-monitor Defibrillator

Model: [CCD-8000C](#)



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26 ,Tianhe,Guangzhou, China





Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguang South Road No 26, Tianhe, Guangzhou, China

Physical characteristics:

Main unit: Portable

Dimensions: L320xW205xH410

Weight: 9.8kg

Technical specifications

DEFIBRILLATOR.

(1) Defibrillator type Manual, Synchronized, Asynchronized,

(2) Energy accuracy $\angle +1\%$

(3) Energy selection 0, 3, 5, 7, 10, 20, 30, 50, 100, 200, 300, 360 joules (nominal at a resistance of 50)

(4) Charging time Maximum 10 seconds to 360 Joules

(5) Standard Adult/Pediatric Paddles Reusable external adult paddles (pediatric paddles integrated)

.Power requirements.

(a) AC Power supply 100V~240V, 50Hz/60 Hz

(b) Vehicle Voltage (Option) DC 12V

Internal Battery Backup

(1) Battery type Rechargeable lead-acid battery (12V 2.0Ah X2)

(2) Charging time Time-Minimum of 4 hours for full charging

(3) Battery charging During AC Power Supply operation of the device.

(4) Charging /Low Battery Indication Red LED

(5) Battery Capacity 100% of the battery will provide 120-min monitoring; and 35 discharges at 360 joules.

Display

-LCD Dimensions 7-inch

-Type High resolution true-color graphic

-Resolution 480X234

ECG

ECG lead: I, II, III, AVR, AVL, AVF, V1~V6

S-T segment analysis: -2.0~2.0mv

Arrhythmia analysis: yes

Protection: With stand 4000vAc/50Hz voltage in isolation and work against electrosurgical knife and defibrillation

Respiration

Measurement method: Thoracic impedance

NIBP

Working mode: Manual/Automatic/stat(5 minutes)

Method: Automatic oscillation

Cuff type: standard Adult; optional Pediatric/Neonatal

Temperature



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguang South Road No 26, Tianhe, Guangzhou, China

Measurement range 20~45°C

Resolution 0.1°C

Pulse Oximetry(spO₂)

Display : spo₂ value, pulse histogram, waveform, Pulse

Range: 0~99%for adult, pediatric and neonatal patients

Probe : standard adult finger clip

optional pediatric Y-type clip and neonatal wrap

Pulse range

Pulse range: 0-255bpm

Standard parameter & configuration

ECG/ RESP/SPO₂/TEMP/PR/NIBP/DEFIBRILLATOR

Optional parameter

Recorder, Li-ion battery (11.1V 4Ah X2)