



Description:

Type: Automatic External Defibrillator

Self-test: Daily, Weekly, Monthly

Mode: Adult, Child

Waveform: Biphasic truncated exponential

Energy: 200 Joule max.

Energy Sequence: Programmable : 30,30,50Joule/ 50,50,75 Joule/120,120,150 Joule /150, 150, 200 Joule

Charge time:

(New, at 25 °C) Less than 6 sec. to 150J

Less than 8 sec. to 200J

Voice Prompt: Extensive voice prompt

Visual Indicators: LED prompts

Control: Two buttons: ON/OFF, Shock

ECG storage: 1500 events.

Data transmission Infrared

Battery

Power: 12V, 3000mAh

Type: Non-rechargeable Li-MnO₂ Cell

Packing:

Dimensions Depth:256mm, Width: 220mm ,Height:(L:65mm, H:80mm)

Weight (with battery): 1.9Kg (4.2 Pounds)