



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

Biphasic Defi-monitor Defibrillator

Model: [CCD-8000A](#)



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Dagan South Road No 26 ,Tianhe,Guangzhou, China





Skype: chinacaremed

Whatsapp: +86 139 2648 0041

Facebook: China care medical

Email: joseph@chinacaremedical.com

Wechat: joseph20130708

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Dagan South Road No 26, Tianhe, Guangzhou, China

Non-synchronizer: outlife difibrillator

Energy: 0、20、50、100、160、250、300、360 joule

Charging time : At 360J less than 10 sec

Paddle options : Reusable external adult and pediatric paddles

Operates from the mains (AC) and with its internal Lead-acid battery

Physical characteristics:

Main unit: Portable

Dimensions: L320xW205xH410

Weight:9kg

Technical specifications

DEFIBRILLATOR.

(1)Defibrillator type Manual, Asynchronized

(2)Energy accuracy $\leq +1\%$

(3) Energy selection 0, 3, 5, 7, 10, 20, 30, 50, 100, 200, 300, 360 joules (nominal at a resistance of 50)

(4)Charging time Maximum 10 seconds to 360 Joules

(5) Standard Adult/Pediatric Paddles Reusable external adult paddles (pediatric paddles integrated)

Power requirements

(a) AC Power supply 100V~240V,50Hz/60 Hz

(b) Vehicle Voltage (Option) DC 12V

Optional parameter

Li-ion battery (11.1V 4Ah X2)