



**Skype:** chinacaremed

**Facebook:** China care medical

**Wechat:** joseph20130708

**Whatsapp:** +86 139 2648 0041

**Email:** joseph@chinacaremedical.com

**Website:** www.chinacaremedical.com

**Add:** Room B205, Changsheng Business Building, Dagan South Road No 26 ,Tianhe,Guangzhou, China

## Defibrillator Trainer AED Machine with CE

### CCS-AED-T



### Specifications:

The CCS-AED-T Trainer is a training tool for preparing emergency responders to use the AED7000 defibrillator to treat a victim of sudden cardiac arrest. The Trainer provides realistic simulation of the defibrillator without actually charging or discharging electrical energy. Voice prompts, simulated energy delivery, pauses for CPR, Training Pads, and the functional system status indicator all accurately mimic the operation of the defibrillator and provide a realistic demonstration.

1. Nine different scenarios are available.
2. Multiple language selections are available.
3. Hand-held IR remote control allows instructors to change training sessions.

### Technology Parameters:

Batteries (Trainer)	CR123A-4x2 12V d.c. 2.8Ah
Battery life (Trainer)	Approximately 200 hours
Battery (remote control)	3V Lithium (CR2025)
Battery life (remote control)	Approximately two years
Number of training scenarios	9
Languages	Chinese English Turkish
Storage temperature without batteries	-30°C to 70°C
Operating temperature	0°C to 50°C
Weight without batteries	1.0 kg
Weight with batteries	1.6 kg



**Skype:** chinacaremed

**Whatsapp:** +86 139 2648 0041

**Facebook:** China care medical

**Email:** joseph@chinacaremedical.com

**Wechat:** joseph20130708

**Website:** www.chinacaremedical.com

**Add:** Room B205, Changsheng Business Building, Daguan South Road No 26 ,Tianhe,Guangzhou, China

---

Dimensions

11.93 X 8.50 X 3.50 in (303 x 216 x 89 mm )