

Dental Soft Tissue Laser



7W/980nm - 10,000Hz

Laser Type	Dental Diode Laser
Model	DD-07
Max Power (CW)	7W
Wavelength	980nm
Operation Mode	CW or Pulsed
Pulse Duration	50µs-1s
Repetition Rate	10,000Hz
Delivery System	200µm, 300µm, 400µm, 600µm fibers
Pilot Beam	3MW/635nm Red Diode
Voltage	AC 110V~220V±20% 50Hz~60Hz±1Hz
Control Mode	5" Full Touch Screen
Activation	Wireless Foot Pedal
Dimensions	155*93*35mm
Weight	0.8kg

