

Rebound Tonometer

Specification:

Range : 3mmHg ~ 70mmHg

Deviation:

$n \pm 1.5 \text{mmHg} (3 \text{mmHg} \leq \text{Pressure} \leq 25 \text{mmHg})$;

$n \pm 2.5 \text{mmHg} (25 \text{mmHg} < \text{Pressure} < 70 \text{mmHg})$;

Characteristic:

➤ Convenient



➤ No Anesthetic



➤ Vertical Horizontal mode



➤ Wireless translate



➤ One-time Probe



Rebound Tonometer

