



LINK CT-600

Chair and Stand Unit

New Technologies And Powerful Functions

Simple operation interface provides
with accurate operation experience.



LINK CT-600

Chair and Stand Unit

SHANGHAI LINK INSTRUMENTS



Features

- ◆ This machine uses American style, it is very smart and fashionable.
- ◆ Table top is easy to rotate and to go up and down.
- ◆ CT-600 is an ideal table unit to install Auto Chart Projector, Slit lamp, Phoropter, Ophthalmoscope and Retinoscope.
- ◆ This table is suitable for hospitals, Ophthalmic Clinics and Optical stores

Specification

Unit size: 1200×800×2100mm

Table Size: 400×500mm

Swing arm oscillation angle: $\pm 35^\circ$

Table up and down range: 810~970mm

Lamp: 220V, 25W

Electric Voltage: AC220V/110V 50/60HZ

Weight: 110kg

Better Detail For Better Quality
Better Detail For Better Quality



Tel: +86-21-69512871
Fax: +86/21-69512872
E-mail: sales@linkcsh.com
Http://www.shlink.en.alibaba.com

SHANGHAI LINK INSTRUMENTS CO., LTD