



# Chair and Stand Unit

New Technologies And Powerful Functions

Simple operation interface provides  
with accurate operation experience.

# Chair and Stand Unit



## Features

- ◆ This machine uses American style, it is very smart and fashionable.
- ◆ Table top is easy to rotate and to go up and down
- ◆ is an ideal table unit to install Auto Chart Projector, Slit lamp, Phoropter, Ophthalmoscope and Retinoscope.
- ◆ This table is suitable for hospitals, Ophthalmic Clinics and Optical stores

## Specification

---

Unit size: 1200x800x2100mm

---

Table Size: 400x500mm

---

Swing arm oscillation angle:  $\pm 35^\circ$

---

Table up and down range: 810~970mm

---

Lamp: 220V, 25W

---

Electric Voltage: AC220V/110V 50/60HZ

---

Weight: 110kg

Better Detail For Better Quality

BETTER DETAIL FOR BETTER QUALITY

