

**Better Detail For Better Quality**  
Better Detail For Better Quality



## **Slit Lamp** **Microscope**


New Technologies And Powerful Functions

Simple operation interface provides  
with accurate operation experience.

# Slit Lamp Microscope

## 14mm Slit Height

14mm aperture provides broader and clearer field of view.




## Metal Joystick

Joystick is made by metal, and the overall appearance is fantastic, and at the same time it increases clinical operation space.



## Five-step drum magnification

Five-step drum magnification of eyepiece. Sclera vein could be observed under the 40X magnification.



## International standards

All the joints can be connected by accessories with international standards.



## Adjustable photoelectric box

By stepless adjustable photoelectric box, it can provide all kinds of appropriate brightness of the lightening according to the clinical doctors' actual demand.

## Filter selection

Four different filters can be inserted into the illumination pathway. Anti-ultraviolet radiation filter ensures patients' safety during examination.

## Specification

Microscope	
Microscope type	Galilean Stereoscopic Microscope
Objective Magnification	Five-step Drum Magnification
Eyepiece	12.5X
Magnifications	6.3X,10X,16X 25X,40X
Diopter Adjustment	-5D~+5D
Pupil Distance	55mm~75mm
Field of view	6.3X=33mm 10X=22.5mm 16X=14mm 25X=8.8mm 40X=5.5mm
Illumination Assembly	
Slit Width	0mm~14mm adjustable
Slit Height	1mm~14mm adjustable
Slit Angle	0°~180° adjustable
Slit Inclination	5°,10°,15°,20°
Diameter of Light Sport	φ0.2mm,φ1mm φ2mm,φ5mm φ10mm,φ14mm
Filters	Heat Absorption,Grey, Redfree,Cobalt Blue
Illumination Source	12V/30W Halogen Source
Power Supply	
Input	AC 110/220V 50Hz/60Hz

