




Rebound Tonometer


New Technologies And Powerful Functions
Simple operation interface provides
with accurate operation experience.

Rebound Tonometer



Specification

Power Supply: DC 3.0V
Battery: Two AAA batteries (batteries not included for export)
Measurement Range: 3 mmHg – 70 mmHg
Measurement Resolution (Tonometry): 1 mmHg
Infrared Transmission Distance: Up to 2 meters
Bluetooth Printer: Thermal printer
Paper Roll Diameter: ≤ φ33 mm
Base: None
Gross Weight: 3.1 kg
Dimensions: 280 × 220 × 110 mm



Better Detail For Better Quality
 Better Detail For Better Quality