



Slit Lamp

New Technologies And Powerful Functions
Simple operation interface provides
with accurate operation experience.

·LIGHT SOURCE SYSTEM UPGRADE

Warm LED light source, suitable for eye tissue for 20,000 hours

Use PWM to ensure the brightness of light source stable

Continuous illumination control makes it more convenient to operate



Features

- PARALLEL TYPE, GALILEAN SYSTEM

Stereoscopic binocular fusion

Can upgrade to teaching slit lamp and digital slit lamp

- MAGNIFICATION CHANGER: REVOLVING DRUM

SL-520: 6x, 10x, 16x, 25x, 40x, 5-step magnification.



Specification

ITEM	SL520
Microscope Type:	Galilean System
Magnification Changer:	Revolving Drum
Eyepiece Magnification:	12.5X
Slit Width:	0mm~14mm Continuously Adjustable
Slit Length:	1mm~14mm Continuously Adjustable
Slit Angle:	0~180°Continuously Adjustable
Filter:	Heat Absorption, Grey, Redfree (Green), Cobalt Blue
Light Control:	PWM: Pulse Width Modulation
Light Source:	Warm LED
Magnification:	6.3X,10X,16X,25X,40X
Diopter:	±7D
Pupil Range:	52mm–85mm
Visual Field Diameter:	6X:33mm; 10X:22mm; 16X:14mm; 25X:8.5mm; 40X:5.5mm
Illumination Type:	Haag–Streit Type
Illumination Tilting:	5°,10°,15°,20°
Slit Aperture Size:	φ0.2, φ1, φ3, φ5, φ10, φ14 (mm)