

Multi-purpose Delivery Bed



Model	Length	Width	Tabletop lowest height	Elevation range
	$\geq 1950\text{mm}$	600mm	$\leq 650\text{mm}$	250mm

The table is hydraulically lifted, The leg plate is removable, various colors of mattress available. The base and column cover all made by premium stainless steel. Good performance on anti-rust and easy to clean.

Parameter

Movements		Model: CCOTYFC2003A
Tabletop trendelenburg	Reverse Trendelenburg	$\geq 10^\circ$
	Trendelenburg	$\geq 25^\circ$
Back plate	Up	$\geq 60^\circ$
	Down	$\geq 10^\circ$

Disposition

Anaesthesia screen	1 piece	Shoulder support (cushion)	2 piece	Arm support (cushion)	2 piece
Handle	2 piece	Leg support (cushion)	2 piece	Foot Support (cushion)	2 piece
Long Fixing Clamp	2 piece	Fixing Clamp	7 piece	Mattress	1 piece

Packaging and delivery

Packaging	Wooden case
Size	142*82*96cm
G./N.W.	195/150 kg