

Electric Operating Table



型号 Model	台面长度 Length	台面宽度 Width	台面最低高度 Tabletop lowest height	台面升降行程 Elevation range
	≥2050mm	550mm	≤710mm	250mm

电动前后倾斜，左右倾斜，背板上下，台面升降，台面平移。底座和立柱罩全部采用优质不锈钢制造。防锈性能好，易清洗。

Table lifting, trendelenburg, reversed trendelenburg, lateral tilt motions, back plate and Horizontal Sliding up folding these motions are controlled by buttons. Electric lifting. The base and column cover all made by premium stainless steel. Good performance on anti-rust and easy to clean.

基本参数 Parameter

动作项目 Movements		型号 Model: YFD04
台面前后倾 Tabletop trendelenburg	前倾 Reverse Trendelenburg	$\geq 25^\circ$
	后倾 Trendelenburg	$\geq 30^\circ$
台面左右倾 Lateral Tilt (Left And Right)		$\geq 20^\circ$
头板上下折 Head plate	上折 Up	$\geq 45^\circ$
	下折 Down	$\geq 90^\circ$
背板上下折 Back plate	上折 Up	$\geq 65^\circ$
	下折 Down	$\geq 30^\circ$
腿板上下折、外折 Leg plate	上折 Up	$\geq 10^\circ$
	下折 Down	$\geq 90^\circ$
	外折 Outward	$\geq 90^\circ$
台面纵向水平移动 Horizontal Sliding Range		$\geq 350\text{mm}$
腰桥上升高度 Kidney Bridge Elevation		$\geq 120\text{mm}$

标准配置 Disposition

麻醉屏架 Anaesthesia screen	1 piece	支肩架 Shoulder support (cushion)	2 piece	搁手架 Arm support (cushion)	2 piece
支身架 Body support (cushion)	2 piece	托腿架 Leg support (cushion)	2 piece	足板 Foot plate (cushion)	2 piece
长滑块 Long Fixing Clamp	2 piece	滑块 Fixing Clamp	9 piece	床垫 Mattress	1 piece

包装数据 Packaging and delivery

包装 Packaging	木箱 Wooden case
尺寸 Size	144*77*101 cm
毛/净重 G./N.W.	210/170 kg