

Electric Gynecology delivery bed



- 1、外围尺寸 1450*600*580-1000mm
Overall size 1450*600*580-1000mm
- 2、高度调节 580-1000mm
Height adjustment 580-1000mm
- 3、背板上折+70°/下折-20°
Back board adjustment degree up +70° /down -20°
- 4、座板上折 0-40°
Seat board adjustment degree up 0-40°
- 5、双面静音轮 75mm，4 个带刹车
Double-sided noiseless wheels 75mm, 4 with brakes
- 6、可伸缩式滑轨带污物盆
Retractable slide with baby basin
- 7、安全承重 180kg
Safe loading weight 180kg
- 8、捷昌电机控制系统
Jiechang Motor Control System
- 9、聚氨酯发泡床面
Polyurethane foam bed surface
- 10、冷钢静电喷涂底座，壁厚 2.0mm
Cold steel electrostatic spraying base, tube thickness 2.0mm