

## Multi-purpose Operating Table



型号 Model	台面长度 Length	台面宽度 Width	台面最低高度 Tabletop lowest height	台面升降行程 Elevation range
	≥2050mm	480mm	≤780mm	250mm

脚踏液压升降。底座和立柱罩全部采用优质不锈钢制造。防锈性能好，易清洗。

Hydraulic lift by foot pedal. The base and column cover all made by premium stainless steel. Good performance on anti-rust and easy to clean.

**基本参数 Parameter**

动作项目 Movements		型号 Model: YFS3001
台面前后倾 Tabletop trendelenburg	前倾 Reverse Trendelenburg	$\geq 30^\circ$
	后倾 Trendelenburg	$\geq 20^\circ$
台面左右倾 Lateral Tilt (Left And Right)		$\geq 20^\circ$
头板上下折 Head plate	上折 Up	$\geq 45^\circ$
	下折 Down	$\geq 90^\circ$
背板上下折 Back plate	上折 Up	$\geq 75^\circ$
	下折 Down	$\geq 20^\circ$
腿板下折 Leg plate		$\geq 90^\circ$
腰桥上升高度 Kidney Bridge Elevation		$\geq 120\text{mm}$

### 标准配置 Disposition

麻醉屏架 Anaesthesia screen	1 piece	支肩架 Shoulder support (cushion)	2 piece	搁手架 Arm support (cushion)	2 piece
支身架 Body support (cushion)	2 piece	托腿架 Leg support (cushion)	2 piece	足板 Foot plate (cushion)	2 piece
长滑块 Long Fixing Clamp	2 piece	滑块 Fixing Clamp	9 piece	床垫 Mattress	1 piece

### 包装数据 Packaging and delivery

包装 Packaging	木箱 Wooden case
尺寸 Size	120*72*104 cm
毛/净重 G./N.W.	185/150 kg