

## Electric Gynecology delivery bed



Type and model	Electric Gynecology Examination Table Obstetric Delivery Bed for Woman Giving Birth
Material:	Stainless steel, power coating treatment is available
Table size	1900×600 mm
Min table height	700 mm
lifting motion	250 mm
Trendelenburg&reversed trendelenburg	≥25°
Back plate up folding	≥70°
Leg plate	Hidden Type
Accessories	Table body/Mattress/Leg support/Arm support/Handle/Anesthetic screen holder/Dirt basin/Remote control