

## Ordinary Operating Table



型号 Model	台面长度 Length	台面宽度 Width	台面最低高度 Tabletop lowest height	台面升降行程 Elevation range
	$\geq 2000\text{mm}$	480mm	$\leq 720\text{mm}$	250mm

各种位置都可以快速、安全地调整。具有经济、实用的特点。底座和立柱罩全部采用优质不锈钢制造。防锈性能好，易清洗。

Various positions can be quickly and safely adjusted. It is characterized by economy and practicability. The base and column cover all made by premium stainless steel. Good performance on anti-rust and easy to clean.

### 基本参数 Parameter

动作项目 Movements		型号 Model: YFS3002A
头板上下折 Head plate	上折 Up	$\geq 45^\circ$
	下折 Down	$\geq 90^\circ$
背板上下折 Back plate	上折 Up	$\geq 75^\circ$
	下折 Down	$\geq 20^\circ$
腿板下折 Leg plate		$\geq 90^\circ$
腰桥上升高度 Kidney Bridge Elevation		$\geq 120\text{mm}$

### 标准配置 Disposition

麻醉屏架 Anaesthesia screen	1 piece	支肩架 Shoulder support (cushion)	2 piece	搁手架 Arm support (cushion)	2 piece
支身架 Body support (cushion)	2 piece	托腿架 Leg support (cushion)	2 piece	床垫 Mattress	1 piece
长滑块 Long Fixing Clamp	2 piece	滑块 Fixing Clamp	7 piece		

### 包装数据 Packaging and delivery

包装 Packaging	木箱 Wooden case
尺寸 Size	118*72*95 cm
毛/净重 G./N.W.	118/90 kg