



## **VERTICAL PRESSURE STEAM STERILIZER**

### **Characteristics:**

- 1. Indicator light indicates working state**
- 2. Over-temperature & over-pressure auto-protection device**
- 3. Double scale indication pressure gauge**
- 4. Safe protection of water lacking**
- 5. With the below exhaust structure,which is convenient to exhaust the cold air**

**6. Automatically shut off with beep reminding after sterilization**

**7. Fully stainless steel structure**

**8. Easy to operate, safe and reliable**

**9. With two stainless steel sterilizing baskets**

### **Technical Data:**

<b>MODEL</b>	<b>150L</b>
<b>Technical Data</b>	
Chamber volume	150L( $\phi$ 510×740)mm
Working pressure	0.22MPa
Working temperature	134℃
Max. working pressure	0.23 Mpa
Heat average	$\leq\pm 1^{\circ}\text{C}$
Timer scope	0~60min
Adjustment of temperature	105~134℃
Power	6KW /AC220V.50Hz
Overall dimension	670×690×1130(mm)
Transport dimension	730×730×1230(mm)
G.W/N.W	130Kg/111Kg

