

Desktop Steam Sterilizer

Characterstics:

- 1.Indicator light indicates working state.
- 2.For 4~6 minutes rapidly sterilizing.
- 3.Sterilizing temperature and time can be preset.
- 4.Steam-water inner circulation system:no steam discharge, and the environment for sterilizing will be clean and dry.
- 5.Over temperature & over pressure auto-protection.
- 6.Safe protection of water lacking.
- 7.Exhaust the cool air automatically.
- 8.Door safety lock system.
- 9.With three stainless steel sterilizing plates.
- 10.The chamber of the sterilizer is made of stainless steel.
- 11.Automatically shut off with beep reminding after sterilization.



Technical data of the Table Top Autoclave Sterilizer:

MODEL TECHNICAL DATA		
Sterilizing chamber volume	35L (ϕ 300×500 mm)	50L (ϕ 340×550 mm)
Working pressure	0.22Mpa	
Working temperature	134℃	
Adjustment of temperature	105-134℃	
Timer	0-60 min	
Heat average	$\leq \pm 1^\circ\text{C}$	
Power	2.7 KW / AC220V 50Hz	
Sterilizing plate	400×200×30 mm (4 pieces)	500×250×30 mm (4 pieces)
Dimension	730×550×510 mm	770×600×540 mm
Package dimension	820×680×650 mm	880×730×700 mm
G. W/N.W	95/80 kg	112/93 kg