



PORTABLE PRESSURE STEAM STERILIZER

Characteristics:

1. Fully stainless steel structure.
2. Hand wheel type of quick-open door structure.
3. Indicator light indicates working state.
4. Double scale indication pressure gauge.
5. Time scope: 0~60min.
6. Automatically shut off with beep reminding after sterilization.
7. Electric heated.
8. Easy to operate, safe and reliable.

Technical data:

Technical Data		
sterilizing volume	18L (280*260mm)	24L (280*390mm)
working pressure	0.14-0.16MPa	
working temperature	126°C	
temperature adjust range	105°C-126°C	
timer range	0-60min	
Max. safety pressure	0.165Mpa	
heat average	$\leq \pm 1^\circ\text{C}$	
power	AC220V. 50Hz / 2KW	
dimension	410×410×520 mm	410×410×630 mm
G. W / N. W	18/16 KG	20/18 KG