



## **PORTABLE PRESSURE STEAM STERILIZER**

**(Automatic microcomputer type)**

### **Characteristics:**

1. Fully stainless steel structure.
2. Hand wheel type of quick-open door structure.
3. Digital display of working status,touch type key.
4. Automatically shut off with beep reminding after sterilization.
5. Electric heated.
6. Easy to operate,safe and reliable.

**Technical data:**

<b>Technical Data</b>		
Sterilizing volume	18L (280*260mm)	24L(280*390mm)
Working pressure	0.14-0.16MPa	
Working temperature	126℃	
Temperature adjust range	105-126℃	
Timer range	0-99 min	
Max. safety pressure	0.165Mpa	
Heat average	$\leq \pm 1^{\circ}\text{C}$	
Power	AC220V. 50Hz / 2KW	
Dimension	410×410×630 mm	410×410×750 mm
G.W / N. W	14/12 KG	15/14 KG