

Desktop Steam Sterilizer



Characteristics:

1. Indicator light indicates working state.
2. For 4~6 minutes rapidly sterilizing.
3. Sterilizing temperature and time can be preset.
4. Steam-water inner circulation system: no steam discharge, and the environment for sterilizing will be clean and dry.
5. Over temperature & over pressure auto-protection.
6. Safe protection of water lacking.
7. Exhaust the cool air automatically.
8. Door safety lock system.
9. With three stainless steel sterilizing plates.
10. The chamber of the sterilizer is made of stainless steel.
11. Automatically shut off with beep reminding after sterilization.

Technical data of the Table Top Autoclave Sterilizer:

MODEL TECHNICAL DATA		
Sterilizing chamber volume	35L (ϕ 300×500 mm)	50L (ϕ 340×550 mm)
Working pressure	0.22Mpa	
Working temperature	134℃	
Adjustment of temperature	105-134℃	
Timer	0-60 min	
Heat average	$\leq \pm 1^\circ\text{C}$	
Power	2.7 KW / AC220V 50Hz	
Sterilizing plate	400×200×30 mm (4 pieces)	500×250×30 mm (4 pieces)
Dimension	730×550×510 mm	770×600×540 mm
Package dimension	820×680×650 mm	880×730×700 mm
G. W/N.W	95/80 kg	112/93 kg