



Technical Parameters

	Model	
Optical	Optics	Achromatic optical system
	Eyepieces	10X / 12.5X Wide-angle lens (optional)
	Binoculars	30°Fixed binoculars, 0°~180°Foldable (optional)
	Pupillary distance	55~75mm
	Diopter	+5D~-8D
	Zoom system	4X~20X Continuous zoom (electric)
	Operating distance	200mm
	Focus	40mm
	Inclination	0~30°
	Assistant Microscope	Modality
Magnification		4X, 7X
Position		Interchange (left/right)
X-Y Focus Function	Movement range	50mm × 50mm
	Reset	Auto-reset
	Speed	2.2mm / sec
Illuminating System	Modality	Aspherical design
	Illuminant	LED
	Control mode	Continuou dimming, Panel control
	Illumination	≥80, 000Lux
	Filter light	Macular protection, Green filter, Cobalt blue (optional)
	Red reflex device	Red color boost
Machine Case	Arm extension length	1250mm
	Up and down adjust	±260mm
	Base size	560mm × 560mm
Multifunctional	X-Y Focus function, Electronic focus, Electronic zoom,	
Foot Control	Illumination adjustment, Illumination switch	
Panel	HD Video system, 4K Video system, 3D Stereo image system	
Other Accessories	Stereo indirect system\Stereo optical invertor\Laser protective goggles	

NO. 1120240101
 With the development of technology, if there are design and parameter changes, please wait for the relevant personnel to inform.



---SURGICAL MICROSCOPE---

OPHTHALMOLOGY


Excellent Optics from OLYMPUS

China Care Medical Equipment Co.,Ltd
 Adress:Room B101, Zhongbang Business Building, PuXing
 Road No. 87, Tianhe, Guangzhou,
 Chinaweb:www.chinacaremedical.com

A stereo & colorful view from Olympus® optics and comfort view eyepiece.

Key feature

OLYMPUS Higher resolution
Excellent depth of view



Independent zoom adjustable focus
(Stereoscopic Same Light Path)



Optics



0°~180°Foldable (optional)




Stereo optical inverter

Stereo indirect system
High resolution
Wide-angle (view)
Undistorted surgical view

Aspherical condenser system

Foot control panel



Red reflex device
Bright surgical view



Video system

Optional video clarity



2D&3D Video solutions
Optional photo and video interface upgraded to full configuration with imaging unit.

Provide better products and services to you with our professional and high-tech