

## Model Tissue Baking Processor



Model:

Voltage: 220V, 50Hz

Power: 300VA

Ambient temperature: (5-40) °C

Relative humidity:  $\leq 85\%$

Range of setting temperature: Ambient  $-80^{\circ}\text{C}$

Temperature Precise error:  $\leq \pm 2\%$

Heating time to reach setting temperature:  $\leq 40$  minutes

Dimension overall: 350\*305\*105mm

External tank dimension: 230\*180\*53mm

N. W. : 3.5kg