

Arm Blood Pressure Monitor



China Care Medical Equipment Co.,Ltd
Address:Room B101, Zhongbang Business
Building, PuXing Road No. 87, Tianhe,
Guangzhou,
Chinaweb:www.chinacaremedical.com

Arm Blood Pressure Monitor

PRODUCT SPECIFICATIONS

External Dimensions (Approximate)	Width: 4.21 in (107 mm) Length: 4.69 in (119 mm) Height: 2.05 in (52 mm)
Weight (Approximate)	0.5 lbs (227 g) without batteries
Display	Digital LCD 41 mm x 61 mm, pressure and pulse displayed
Power Options	4 AAA alkaline batteries, 6VDC AC power supply Input: 100-240 V, 50-60 Hz, 350 mA Output: 5 V, 1 A
Method	Oscillometric
Measurement Range	Rated cuff pressure: 0 mmHg - 300 mmHg (0 kPa - 40 kPa) Systolic Measurement Pressure: 55-255 mmHg Diastolic Measurement Pressure: 25-200 mmHg Pulse Measurement Pressure: 40 to 199 beats/minute
Accuracy	Pressure: +/- 0.4kPa (3 mmHg); Pulse Value: +/- 5%

PARTS AND ACCESSORIES

CONFIGURATION

Arm Blood Pressure Device, Size 22-32cm BP cuff **FLEXIPOINT BL**

OOD PRESSURE CUFFS

Soft cuff: Child (17–22 cm)

Soft cuff: large Adult (22–42 cm)

Rigid cuff: large Adult (22–42 cm)

ACCESSORIES

Adapter: Optional US, EU, UK, AU adapters

Voice: Optional Custom Languages

Battery: Optional 4*AAA Alkaline batteries

China Care Medical Equipment Co.,Ltd
Address: Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
China web: www.chinacaremedical.com

