

1. Drag 3



- Use a stainless steel air intake pipe and air intake fitting
- Avoid oxidation and rusting of the joint intake pipe
- Adopt stainless steel intake pipe and intake joint
- Avoid oxidation and rust of joint inlet pipe

Product Parameters

Available	1-3 treatment machines
Using voltage	AC220/50Hz
Rated power	1120W
Speed: r/min	1380 r/min
Flow rate: L/min	140 L/min
Volume L	60 L
Starting pressure Starting pressure (MPa)	0.5Mpa
Exhaust pressure	0.8Mpa
Noise level	70dB (A)
External dimensions: Boundary dimension (mm)	800×350×690
Weight	66kg



Optional: Automatic drainage for worry-free use

Optional: Automatic drainage so that you can use at ease



Motor standalone start-stop function
Motor separate start and stop function



The motor features shock-absorbing springs for quieter operation.
Motor with shock-proof spring, more quiet



Rust-proof treatment inside the tank container
Anti-rust treatment inside tank