



### GENERAL FEATURES

- LCD display screen, 6.0 inches.
- Reset feature automatically resumes all settings back to factory default options.
- IP54 waterproof.
- 1-2 points calibration with Standard Recognition.
- Manual temperature compensation ensures accurate results

Model		
Parameters		pH/mV
pH	Range	0.00 to 14.00 pH
	Resolution	0.01pH
	Accuracy	± 0.01 pH
	Calibration Points	Up to 2
	Standard Recognition	NIST buffers
mV	Range	-1800 to 1800 mV
	Resolution	1 mV
	Accuracy	± 0.1%FS
Measurement	Reading Mode	Continuous
	Temp. Compensation	MTC
Inputs	pH Electrode	BNC(Q9)
Display options	Backlight	Yes
	Auto Shutdown	300, 600, 1200, 1800, 3600 sec., off
	IP Rating	IP54
General	Power	AC Adapter, 100-240V AC input, DC9V output
	Dimensions	200 × 160 × 63 mm
	Weight	600g (1.32 lb)

China Care Medical Equipment Co.,Ltd  
 Address: Room B101, Zhongbang Business Building, PuXing  
 Road No. 87, Tianhe, Guangzhou,  
 Chinaweb: [www.chinacaremedical.com](http://www.chinacaremedical.com)