



Microplate Washer

China Care Medical Equipment Co.,Ltd  
 Adress:Room B101, Zhongbang Business Building,  
 PuXing Road No. 87, Tianhe, Guangzhou,  
 Chinaweb:www.chinacaremedical.com



Microplate Washer

**Feature:**

- One or more strips of ELISA Plate washing key board choice directly
- Minimal residual volume (<1  $\mu$  L/well)
- Dispense precision Cv<2%
- Detachable ELISA Plate pallet, convenient for the pallet clearing and sterilization
- Inter 100 groups wash programs, Adjustable and programmable washing modes
- Cleaning head 5-locations location adjusted,
- Soaking or vibration time programme
- Temperature:10  $^{\circ}$ C -30  $^{\circ}$ C wet  $\leq$ 80%
- Power Supply: 220V110VAC  $\pm$  10%, 50/60Hz 100W
- Automatic waste and lotion and distilled water liquid warning



**Specification:**

Micro-plate Types: 48 / 96-well U-V- / Flat-bottom plates or strips  
 Dispense volume:0-950 $\mu$ l,adjustable  
 Washing Cycles:1-99 optional  
 Washing Mode: 1-12, row and rank, selectable  
 Soaking time:0-999s, adjustable  
 Vibration time :0-999s,adjustable  
 Display: Large LCD  
 Channel: 2 channels (including one channel for water)  
 Memory: Up to 100 programs can be stored  
 Interface: Rs232  
 Power supply:220V  $\pm$  10% 50/60Hz AC110V  $\pm$  10% 50/60Hz 100W  
 Dimension: 440mm (L)  $\times$  305mm (W)  $\times$  154mm (H)  
 Net Weight: 13kg



Microplate Washer

Built-in internal-temperature  
 cyclic ELISA Plate pallet incubator(patent)  
 temperature precision 37  $\pm$  0.3 $^{\circ}$ C

