



330lb Load Bearing



Foldable for easy transport



sturdy carbon steel material



Shock absorption springs



Product parameters

整车尺寸	Overall Dimension	43.3"×24.4"×37" / 110*64*95cm
折叠尺寸	Folding Size	31.8"×14.6"×28.7" / 81*37*73cm
座椅尺寸	Seat Size	17.3"×18.1" / 44*46cm
整车净重	Net Weight	93.7lb / 42.5kg (including 12AH lead-acid battery)
最大负载	Maximum Load	330lb / 150kg
电池	Battery	24v-12/20ah Lead-acid Battery 12ah 8.5kg / 18.7lb 20ah 13kg / 28.7lb
电机	Motor	24v 250w×2
续航里程	Driving Distance	15-25km
爬坡能力	Climbing Capacity	≤9°
最大前进速度	Max Forward Speed	6km
前轮规格	Front Wheel Size	10"
后轮规格	Rear Wheel Size	16"
包装尺寸	Package Size	32.2"×17"×31.3" / 82*43*79.5cm