

Product Characteristics

Homogenizer



◆ Large capacity, high power, high rotational speed and simple operation

The homogenizing cup for crushing is made of stainless steel, which has good corrosion resistance and is convenient for cleaning and disinfection

Common tools for laboratories, analysis rooms, education and scientific research

◆ Experimental equipment applied to the mixing and crushing of liquid media and soft solids

Product name	Tissue Triturator
power	200W
Capacity	1000ml
Speed range	start up ~12000 rpm
Voltage	220V 50Hz
Dimension(mm)	230x230x420
Package(mm)	280x280x470
G.W./N.W.	8/7kg

China Care Medical Equipment Co.,Ltd
 Address:Room B101, Zhongbang Business Building,
 PuXing Road No. 87, Tianhe, Guangzhou,
 Chinaweb:www.chinacaremedical.com