

**Product Characteristics**

**Analysis balance FA series**



- Electromagnetic force compensation technology
- Aluminum alloy die cast, stainless steel plstter
- Super bright LCD display with backlight.
- Clear glass windshield wupplied as standard
- Overload Protection, bubble level adjustment.
- Under hook for huge object.
- Full Capacity subtraction
- Multi weighing unit conversion g/mg/ct/oz
- Stabilization time, 3-5 seconds.
- Check weighing, parts counting, percentage weighing.
- Universal power adaptor supplied as standard

Weighing scope	0-200g
Precision	0.1mg
Repeatability	±0.2mg
Linear error	±0.2mg
Pan size	Φ80mm
Interface inside	RS232C
Dimension	34cm×21.5cm×35cm
Voltage	AC 110-240V
Package(mm)	480*340*500
G.W./N.W.	9/8KG

China Care Medical Equipment Co.,Ltd  
 Adress:Room B101, Zhongbang Business Building, PuXing  
 Road No. 87, Tianhe, Guangzhou,  
 Chinaweb:[www.chinacaremedical.com](http://www.chinacaremedical.com)