



Excellent pressing and dotting mechanism
ensure accurate measuring result for each lens.



SUPORE BRILLIANT

The range of measurement	
CCE-LM-700	
Sphere:	0m-1to±25m-1
Cylinder:	0m-1to±10.00m-1
Astigmatism:	0° to 180°
ADD:	0m-1 to±10.00m-1
Prism basal Angle	0° to 360°
measurement rang	
Prism:	Horizontal 10cm/m, Vertical 15cm/m
The accuracy of measurement	
Dioptric:	(0.01/0.12/0.25)m-1
Prism:	(0.01/0.12/0.25)cm/m
Mode of measurement	
Cylinder:	+, ±, -
Prism:	X-Y, P-B
Speed of measurement:	0.2 second
Diameter od measurable lens:	18 to 110mm
Other parameters	
Weight:	About 5.1kg
Power supply:	AC 220V 50Hz

China Care Medical Equipment Co.,Ltd

Room B101, ZhongBang Business
Building, PuXing Road No. 87, Tianhe,
Guangzhou, China

www.chinacaremedical.com

www.chinacaremedical.com



SUPORE BRILLIANT

AUTO LENS METER



AUTO LENS METER

Concise and compact design
saves installation space.

Tilting screen is convenient
for users to observe.

Intuitive, simple and clear
interface makes operation
more easier.



Product features




Measuring single/binocular and progressive
multi-focal lenses with ease and rapidity.



Simple user operating interface.




LCD liquid-crystal display




Provide more accurate UV measuring system.



Hartmann principle is more accurate for
measurement.



QR Code Printing



Electronic PD Ruler




Strengthening dotting device



It is convenient to replace paper
with integration printer.



Green light measurement without
compensation of Abbe



Automatic R/L eye recognition

China Care Medical Equipment Co.,Ltd
Address:Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
China
web:www.chinacaremedical.com