

Accessory Packing List

Item No.	Description	Quantity	Unit	Remark
203050000014314	First-stage combustion tube	2	pieces	
203050000014315	Secondary combustion tube	3	pieces	
203010000005333	Drying tube	4	pieces	
f01000000018965	Linear copper oxide	6	bottles	
f01000000018967	Linear copper	2	bottles	
f01000000018963	Tungsten grain bracket	1	box	
f01000000018969	Filling ball	3	bottles	Ball diameter 3-5mm
f01000000018966	Platinum catalyst	4	bottles	
f01000000018968	Silver	1	pack	
f01000000018970	Desiccant	1	bottle	
f01000000018972	Tin foil	5	boxes	
f01000000018964	Aspartic acid	1	bottle	
f01000000018971	Filter fillers	1	pack	
103050000000001	Power cord	1	piece	
104010000005347	Serial line	1	piece	
104010000000051	USB re-serial line	1	piece	
204060000014569	Fluoroelastomer seal (inside diameter x wire diameter).	15	pieces	φ9x3
203040000005338	Silicone seal	8	piece	φ16x3.1
203040000005339	Sampler seal	2	pieces	
203040000006984	Fluoroelastomer seal	6	pieces	φ24x1.5mm
203040000006985	Fluoroelastomer seal	6	pieces	φ27x1.5mm
203040000006980	Silicone seal	6	pieces	φ23x3.1mm
203040000012534	Fluorine seal	15	pieces	φ20.35x1.78
203040000012533	Fluorine seal	6	pieces	φ30.2x2.4
d03120000005362	Grease	1	box	
d01120000009263	Stainless steel scraping	1	piece	
205000000015035	Supplies filling gaugies	1	set	
203050000005354	Restore the pipe access tool	1	piece	
205000000007361	Sample press plate	1	cover	

d02020000005343	High temperature gloves	1	pair	
d01030000005358	Pliers	1	piece	
205000000007357	Gray powder cleaning tool	1	set	
203020000007356	Consumables replacement bracket	1	piece	
204040000007506	Copper pipe	30	meters	
d01060000009223	VETUS tweezers ST-13 flat circle	1	piece	
d01030000009225	Allen key	1	piece	
d01020000011758	Allen key	1	piece	
804040000021172	Accessory packing list	1	piece	
203050000011642	Oxygen gun tools	1	set	
205000000005351	Grey powder tube tongs	1	set	
204000800001338	Oxygen gun set	10	pieces	
203040000007036	Silicone tubes ϕ 2.5x ϕ 1.5	30	cm	
203050000014566	Gray powder tube	2	pieces	
203090000015030	Ceramic crucibles	10	pieces	
d01030000011649	T-type hex wrench	1	piece	
203030000005360	Sample box	1	piece	
203030000005361	Funnel	1	piece	