

Product Characteristics



- ◆ Electromagnetic force sensor
 - ◆ Super bright detachable displayer
 - ◆ Multi-unit weighing function, customizable units
 - ◆ RS232 and USB bidirectional interface
 - ◆ Glass windshield-antistatic material
 - ◆ Fast sensitivity setting
 - ◆ Below weighing system
 - ◆ Time display
 - ◆ The weighing results can be exported to Excel
 - ◆ Weighing function
- Basic weighing, Dynamic weighing, Piece counting, Percent weighing, Density determination, Σ weighing function, Pipette calibration function

MODEL			
Weighing range	100g	120g	220g
Reading accuracy	0.01mg	0.01mg	0.01mg
Stabilization time	3 S	3 S	3 S
Preheating time	30-60 MIN	30-60 MIN	30-60 MIN
Weighing units	g/oz/ct/mg	g/oz/ct/mg	g/oz/ct/mg
Repeatability	±0.03mg	±0.03mg	±0.03mg
Linearity error	±0.05mg	±0.05mg	±0.05mg
Output interface	RS232 standard / USB standard	RS232 standard / USB standard	RS232 standard / USB standard
Scale plate size	Φ90mm	Φ90mm	Φ90mm
Power supply	AC110V-240V	AC110V-240V	AC110V-240V
package	490x380x485mm	490x380x485mm	490x380x485mm
G.W./N.W.	7/6KG	7/6KG	7/6KG

