



Handheld Urine Analyzer



Easy to Operate

One-click Operation



Rapid Test

14 Items in 60s



Intelligent Calibration

Automatic Color Compensation

> Parameters of The Device

Device

Handheld Urine Analyzer BK-150

Detection Principle

Dual Wavelength Reflection Method

Sample Type

Urine


Size

178 mm (length) × 87 mm (width) × 46 mm (height)

Weight

360 g

> Test Item



China Care Medical Equipment Co.,Ltd
Address: Room B101, Zhongbang Business Building, PuXing Road No. 87, Tianhe, Guangzhou, China
web: www.chinacaremedical.com