



Automatic 3 diff Hematology analyzer

Technical parameters

Test items

WBC, LYM #, MID #, GR #, LYM%, MID%, GR%, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-CV, RDW-SD, PLT, MPV, PDW, PCT, P-LCC, P-LCR

Histogram

White blood cell (WBC) histogram, red blood cell (RBC) histogram, platelet (PLT) histogram

Sample type

pulse blood, peripheral blood, prediluted blood

Test mode

whole blood mode, pre-dilution mode

Blood volume

10 µL of venous blood and peripheral blood, 20 µL

Storage

automatic storage, maximum of 200,000 sample results (including sample information)

Throughput

60 samples / hour, can be turned on for 24 hours, custom sleep time

Data query method

Forward and backward flushing and blockage removal technology, combined with high-pressure burning and blockage removal, with extremely low blockage rate

Operating characteristics

Forward and backward flushing and blockage removal technology, combined with high-pressure burning and blockage removal, with extremely low blockage rate

Display

10.4 inch HD touch LCD screen, all results and histograms are displayed on the same screen

Operation mode

support keyboard, mouse operation, graphic touch operation, handwritten input

Quality control

L-J, X-R, floating mean value, etc

Printing

built-in thermal printer, a variety of printing formats, and support for external printer

Information transmission

bidirectional LIS communication, HL7 and other transmission protocols

Wireless connection

WIFI, Bluetooth (optional) and other wireless transmission of data, wireless printing, LIS connection, laboratory LAN

Interface

LAN, RS232, HDMI, 4 USB

Language

Chinese operation, support multi-language

Options

bar code scanner (optional)

Power supply and power

~100-240V

Weight

Net weight: about 17kg Gross weight: about 23kg

Volume

Machine size: 265 * 380 * 420mm
package size: 566 * 433 * 580mm

Distributor:

New upgrade 3 Diff of small volume and big wisdom

Automatic Hematology Analyzer



Automatic 3 diff Hematology analyzer

China Care Medical Equipment Co.,Ltd
Address:Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
Chinaweb:www.chinacaremedical.com

