

玻璃体切割手术器械包

序号	产品单位	数量
1	Intraocular Retinal Micro Forceps, 23G straight membrane forceps with platform tip	把 1
2	Intraocular Retinal Micro Forceps, 23G straight membrane forceps without platform tip	把 1
3	Intraocular Retinal Micro Forceps, 23G straight forceps with 1×2 teeth tip	把 1
4	Intraocular Retinal Micro Forceps, 23G extended ILM forceps tip	把 1
5	Intraocular Scissors, 23G horizontal straight cutting tip	把 1
6	Intraocular Scissors, 23G horizontal curved cutting tip	把 1
7	Intraocular Foreign Body Gripper, 23G two-prong perforated tip	把 1
8	Intraocular Retinal Micro Forceps, two-finger interchangeable handle	把 1
9	Intraocular Foreign Body Gripper, single-finger interchangeable handle	把 1
10	Retinal Hook, 23G Type 1	把 1
11	Retinal Hook, 23G Type 2	把 1
12	Intraocular Membrane Spatula, single curved tip, titanium	把 1
13	Corneal Epithelial Flap Dissector, double-ended Type C, titanium	把 1
14	Flute Needle, standard 23G	把 1
15	Microsurgical Instrument Sterilization Tray, 265×162×28mm	只 1

