

鼓膜切开手术器械包

NO.	SPECIFICATION	QUANTITY
1	12.5cm dissecting scissors, str	1 pc
2	Ear Speculum 4, 5.5, 7-3x1	1 set
3	Cotton Carrier , ϕ 1	1 pc
4	Curette, Ear, 16cm, Straight Round Head ϕ 1.0	1 pc
5	Curette, Ear, 16cm, Corner Bend, 25° Oblong Head, Width 1.0/1.2	1 pc
6	Cushing Bayonet Forceps, 16cm, Serrated	1 pc
7	Knife, Myringotome, 19cm, Triangular Blade, Bayonet Shaped	1 pc
8	Fraser Suction Tube, 22cm, Curved, Φ 2.5mm	1 pc
9	Fraser Suction Tube, 27cm, Straight, Φ 1.5mm, Simple	1 pc
	Casing	1 ea

