


NO.	SPECIFICATION	数量	
1	kidney tray 204×130×25	1	
2	Scalpel handle 4#	2	
3	Scalpel handle 7#	1	
4	14cm tissue forceps, 1x2 teeth	2	
6	20cm tissue forceps, 1x2 teeth	1	
7	Dressing forceps 20cm	1	
8	Scissor, Mayo, Surgical, 16cm, Straight	2	
9	Scissor, Mayo, Surgical, 18cm, Curved	1	
10	Scissor, Metzenbaum, 20cm, Curved	1	
11	Elevator, Periosteum, 28cm, Straight, Round Edge, Heavy Duty Blade, Width 20mm	2	
12	Laminectomy Punch	2	
13	Nucleus Pulposus Forceps(Single Ring Handle) 3X10	1	
14	Retractor, Cervical, 16cm, Movable Single Joint, With 5 Hooks	1	
15	Retractor, Cervical, 12cm, Rack Type Single Joint	1	
17	Brain, Reflex Hammer, 14cm	1	
18	Mallet, Bone, 22cm, Ball Plane, 70g	1	
19	Fraser Suction Tube, 22cm, Straight, Φ5mm, Simple	1	
20	Fraser Suction Tube, 22cm, Curved, Φ3mm	1	

21	Blade, Osteotomes, 23cm, Straight Edge Blade, Width 8, Corner Handle	1	
22	Lambert Rongeurs 1	1	
23	Scissors, Bone, 18cm, Double Joint, Straight Tip	1	
24	Rongeurs, Bone, 22cm, Double Joint, Angular Curved 20° Blade, 2, Non-Angular Shank	1	
25	Rongeurs, Bone, 22cm, Double Joint, Angular Curved 20° Blade, 6, Non-Angular Shank	1	
26	Rongeurs, Bone, 22cm, Double Joint, Angular Curved 20° Blade, 8, Non-Angular Shank	1	
27	Forceps, Wire Cutting 18cm, Single Joint, Tiger Head	1	
28	14cm Backhaus forceps,	4	
29	14cm forceps, Crile, str	2	
30	12.5cm mosquito forceps, cvd	4	
31	14cm forceps, Crile, cvd	4	
32	18cm forceps, Crile, cvd	4	

33	16cm Allis forceps	6	
34	18cm forceps, Kocher, str	2	
35	18cm Needle holder, Mayo-Hegar	2	
36	20cm Needle holder, Mayo-Hegar	2	
37	25cm Sponge holding forceps, str, grooved 8mm	2	

