



China Care Medical Equipment Co.,Ltd
 Address:Room B101, Zhongbang Business Building, PuXing
 Road No. 87, Tianhe, Guangzhou,
 Chinaweb:www.chinacaremedical.com

Trolley Color Doppler

Feature:

- Built in base on windows System
- Hard disk 120GB, support enlargement of memory
- Support many kinds of printers
- Weight: 50 Kg
- Package size:900 x690 x1160 mm
- 15inch LCD medical monitor
- Probe interface:4
- Other types port:2USB
- port,VGA,W-Lan,Video,power
- Adapter,Power cable

Mode:

- B&W mode:B,2B,B/M,4B,M,PW
- Color mode:CFM,PDI,PW,B/C,
- Frequency range:2.0-10.0Mhz

Standard configuration:

- 3.5Mhz Convex probe
- 7.5Mhz Linear probe

Optional:

- 6.5Mhz Vaginal probe
- 5.0Mhz Mirco convex probe



Ultrasound machine

Feature:

- 10inch LCD monitor
- Light weight and thin
- Color:white&blue
- Probe interface:2

Standard configuration:

- 3.5Mhz Convex probe

Optional:

- 6.5Mhz Vaginal probe
- 7.5Mhz Linear probe
- 5.0Mhz Mirco convex probe

