

-25°C Biomedical Freezer



Applications

The -25°C medical deep freezers are designed to store vaccines, blood plasma, reagents, and many other biological materials which are widely used in research institutions, clinics, and medical fields.

Features

- Microprocessor control, ensure good inner temperature uniformity, accuracy is 0.1°C
- HC green refrigerant is energy-saving and environmentally friendly.
- Wide available voltage range: 220V±10%.

Security

- Safety door lock and controller password to prevent unauthorized access.
- Built-in backup battery to display temperature and alarms when power failure or system error.
- Visual and audible alarm system.
- Power failure protection: turn-on delay of the cooling system after power failure.
- Controller error protection: the cooling system will remember the normal working cycle and keep working depending on this memory when the controller or two sensors failed.

Humanized design

- Drawers for easy storage.



Specifications

Main Specifications	
Temperature Range	-10 ~ -25 °C
Ambient Temperature	10 ~ 32°C
Capacity (L/Cu. Ft.)	100 / 3.53
Interior Dimensions (W×D×H) mm	410×410×670
Exterior Dimensions (W×D×H) mm	540×553×877
Weight	32 kgs
Noise	50 dB
Control and Power	
Controller	Microprocessor
Display	Digital
Sensor	NTC
Voltage and Frequency	220V±10%, 50/60Hz
Power	73 W
Electric Current	0.7 A
Power Consumption (kWh/24h)	0.95
Audible and Visual Alarms	
High / Low Temp. Alarm	Yes
Power Failure Alarm	Yes, 8h
Sensor Error Alarm	Yes
Controller Error Alarm	Yes
Door Ajar Alarm	Yes
Refrigeration	
Refrigeration Type	Direct Cooling
Defrost	Manual
Compressor (Brand/pcs)	Chinese Brand / 1
Refrigerant / Weight	R600a / 37g
Cabinet	
Door (pcs)	Foamed Solid Door / 1
Drawer (pcs)	4
Interior Cabinet Material	HIPS
Exterior Cabinet Material	Coated Steel
Insulation Material	PURF (δ78mm)
Door Lock	Yes
Fixed Feet (pcs)	4
Accessories	
Remote Alarm Port	Yes
USB Data Download Port	Optional
Temperature Printer	Optional
RS485 Port	Optional
Temperature Datalogger	Optional
Shipping	
Package	Carton box
Shipping Weight	43 kgs
Shipping Size (W×D×H)	640×650×1100
Shipping Volume	0.46 CBM

* Due to continuous product upgrades, the above specifications are for reference only. For the latest specifications, please contact us.

China Care Medical Equipment Co.,Ltd

Address: Room B101, Zhongbang Business Building, PuXing Road No. 87, Tianhe, Guangzhou,

China web: www.chinacaremedical.com