

-60°C ULT Freezer



Applications

The -60°C ultra-low temperature freezers are designed for the long-term preservation of viruses, germs, vaccines, biological tissues and organs, special food, medicines, reagents, etc.

Features

- Accurate temperature controller, setting temperature range is -25 ~ -60°C, accuracy is 0.1 °C.
- Suggested ambient temperature: +10~32°C.
- Stainless Steel inner chamber for friendly use and easy cleaning.
- Brand compressor to ensure excellent stability of the cooling system.
- Patented refrigeration technology, high-efficiency mixed gas for the cooling system.

Security

- Safety door lock to prevent unauthorized access.
- Visual and audible alarm system: High/low temperature, sensor failure, controller error.
- Power failure protection: turn-on delay of the cooling system after power failure.
- Controller error protection: the cooling system will remember the normal working cycle and keep working depending on this memory when the controller or two sensors failed.

Humanized design

- Casters with stoppers for easy moving and fixation.
- Special-designed service-free drainage system.



Specifications

Main Specifications	
Temperature Range	-25~-60℃
Ambient Temperature	10~32 ℃
Capacity (L/Cu. Ft.)	485 / 17.13
Interior Dimensions (W×D×H) mm	1455×525×665
Exterior Dimensions (W×D×H) mm	1670×810×930
Weight	124 kgs
Noise	53 dB
Control and Power	
Controller	Microprocessor
Display	Digital
Sensor	PTC
Voltage and Frequency	220V±10%, 50/60Hz
Power	600 W
Electric Current	3.44 A
Power Consumption (kWh/24h)	7.45
Audible and Visual Alarms	
High / Low Temp. Alarm	Yes
Sensor Error Alarm	Yes
Controller Error Alarm	Yes
Refrigeration	
Refrigeration Type	Direct Cooling
Defrost	Manual
Compressor (Brand/pcs)	CUBIGEL / 1
Refrigerant / Weight	R4001 / 320g
Cabinet	
Door (pcs)	Foamed Solid Door / 1
Interior Cabinet Material	Stainless Steel
Exterior Cabinet Material	Coated Steel
Insulation Material	PURF (100mm)
Door Lock	Yes
Test Hole (pcs)	1
Leveling Castors (pcs)	4
Accessories	
Temperature Datalogger	Optional
Protective Gloves	Optional
Shipping	
Package	Carton box with pallet
Shipping Weight	149 kgs
Shipping Size (W×D×H)	1750×880×1040 mm
Shipping Volume	1.6 CBM

* Due to continuous product upgrades, the above specifications are for reference only. For the latest specifications, please contact us.

China Care Medical Equipment Co.,Ltd

Address:Room B101, Zhongbang Business Building, PuXing Road No. 87, Tianhe, Guangzhou,
China web: www.chinacaremedical.com