



Technical Parameters:

Image Storage	64 Frames
Display Depth	≥ 160
Scanning Model	Electronic Linear Array, Electronic Convex Array
LCD Size	9.7 Inch/ 12Inch
Gray Scale	256 Levels
Image Flip	Up/ Down, Left/ Right, Blake/ White
Image Processing	Image Smoothing/ Sharpening, Tissue Harmonic, Histogram, DR, Gamma Correction, Pseudo-Color
Measurement	Freehand, Ellips, Distance for Perimeter, Area Volume, Heart, G.A, EDD, BPD-FW, FL, AC, HC, CRL, AD, GS, LMP
Character Display	Date, Clock, Name, PID, Age, Sex, Hospital Name, Doctor
Notation	Full-screen Character Editor, Body Mark, Position Indication
Report	Abdomen and gynecology, Cardiac, Obstetrics, Urology
Display Mode	B, B+B, B+M, M, 4B
Scanning Angle	Adjustable
Cine Loop	≥ 400 Frames
Output	VGA and PAL video outputs
Power Consumption MAX	100VA
Product Size	289 × 304 × 222mm
Carton Size	395 × 300 × 410mm
N.W./ G.W.	6kg/ 7kg