





Full-Digital Color Doppler



Muti-Department

Application



High grade display

- ◎ High resolution 15" LED
- ◎ 80 elements probe
- ◎ Support infinite amplification under real time scan (B, B+C, 2B, 4B)
- ◎ Color blood flow, PW, PDI

Human function &

Keyboard design

- ◎ Optimized keyboard
- ◎ Backlight keyboard

Data store & software

- ◎ DICOM 3.0
- ◎ Internal working station: support hard disk, search & browse patients data
- ◎ Powerful measuring & calculating software for clinical demands